

# Purification of SARS Hyperimmune Globulins

*The 4th Plasma Product Biotechnology Meeting  
9-12 May 2005*

Bing L. Wong, Ph.D.  
Advantek Biologics Limited  
[blwong@advantekAB.com](mailto:blwong@advantekAB.com)



*We save lives through safe biologics*

# Severe Acute Respiratory Syndrome

## Global Outbreak 2002-2003

<b>Countries</b>	<b>Cases</b>	<b>Death</b>	<b>Fatality</b>
China	5,327	349	7%
Hong Kong	1,755	299	17%
Taiwan	346	37	11%
Canada	251	43	17%
Singapore	238	33	14%
Viet Nam	63	5	8%
All countries	8,096	774	9.6%

From Nov 16, 2002 to Jul 13, 2003.

Source: WHO



*We save lives through safe biologics*

# Current First-line Treatment for SARS

- Ribavirin and Steroids

Retrospective study of 323 patients (Princess Margaret & Wong Tai Sin Hospitals, Hong Kong)

64%, rapid disease progression

21%, intensive care

13%, ventilator support

8%, death

Tsui et al *Emerg Infect Dis* 2003; 9: 1064-9



*We save lives through safe biologics*

# Convalescent Plasma Treatment of SARS

- Zhou *et al*, *Zhonghua Yi Xue Za Zhi*. 2003; 83: 1018-22

a 74-year-old Beijing patient received 50 ml

Outcome: complete recovery

- Soo *et al*, *Clin Microbiol Infect*. 2004; 10: 676-8

28 Hong Kong patients received 200-400 ml


Outcome: no death

shorter disease course

good responses before day-16



# World's First SARS Investigational Drug Human Hyperimmune Globulins



2003



*We save lives through safe biologics*

# Collection of Convalescent Plasma

- At least 28 days after hospital discharge
- Plasmapheresis collection of 200-400 ml/unit, (total 47 units), frozen immediately below -20 °C
- Routine screening negative for HBsAg, HCV antibody, HIV antibody and syphilis, and normal for ALT
- Tested positive for anti-SARS antibodies



# Determination of anti-SARS Titer

- ELISA (commercial kit)
- Immuno-Fluorescence Assay (supported by the Chinese University of Hong Kong)
- Neutralizing Antibody Test (provided by the Chinese National Institute for the Control of Pharmaceutical and Biological Products)



# Anti-SARS Titer of Convalescent Plasma

ELISA Titer	Number of Plasma Units
1 : 2	2
1 : 4	5
1 : 8	12
1 : 16	17
1 : 24	6
1 : 32	3
1 : 48	2


ELISA titer of pooled convalescent plasma, 1 : 12



*We save lives through safe biologics*




# Process: Human Hyperimmune Globulins



# Titer of Human Hyperimmune Globulins

ELISA Titer	1 : 83
IFA Titer	1 : 1,600
Neutralizing Titer	1 : 200



# PCR Detection of SARS Virus RNA

- RNA extracted from the human hyperimmune globulins
- 200 ng RNA used in the quantitative real-time PCR reactions, with proper negative, positive and spiking controls, following *Biochem Biophys Res Commun.* 2003; 312: 1290-6
- Result: Negative for SARS virus RNA




# Characteristics of Human Hyperimmune Globulins

Assay items	IVIG specifications	SARS IgG
pH	3.8 – 4.4	4.1
Purity	≥ 95%	99%
Mono + Dimer	≥ 95%	100%
ACA	≤ 50%	11%
PKA (IU/ml)	≤ 35	Not detected
TnBP residue	≤ 10 ppm	< 2 ppm
Triton residue	≤ 10 ppm	< 5 ppm



# Process for Equine F(ab')<sub>2</sub> Fragments



# Equine Serum with SARS Antibodies

- SARS virus inactivated with formaldehyde, inactivation confirmed in cell culture
- Immunization/re-immunization of horses
- Much higher neutralizing titer 1:10,280, compared with human convalescent plasma
- Negative for SARS virus RNA, with the quantitative real-time PCR



## Characteristics of SARS F(ab')<sub>2</sub>

Assay items	Antitoxin specifications	SARS F(ab') <sub>2</sub>
Protein conc.	≤ 17%	2.1%
F(ab') <sub>2</sub> Purity	≥ 60%	64.3%
pH	6.0 – 7.0	6.59



# Comparison: Human v.s. Equine

## Human SARS Globulins

- Very limited
- Lower titer
- Less antigenic

## Equine SARS F(ab')<sub>2</sub>

- Unlimited
- Much higher titer
- More antigenic





# Conclusion

- Purified and concentrated SARS hyperimmune globulins are prepared and ready to be further evaluated.
- The SARS hyperimmune globulins may help to better contain future SARS outbreaks



# Acknowledgment

Shenzhen Weiwu Guangming Biological Products Co Ltd

Mr. Jiling Hong, Mr. Xin Zhang, Mr. Zhan Zhang

The Chinese University of Hong Kong

Prof. John S. Tam, Prof. Gregory Cheng

Sun Yat-sen University

Prof. Jiahai Lu



*We save lives through safe biologics*